CAMP[®]
SCANDINAVIA

AFLASTNING

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OPTIMA[®]
molliter

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Emerging Treatments and Technologies
ORIGINAL ARTICLE

An Off-the-Shelf Instant Contact Casting Device for the Management of Diabetic Foot Ulcers

A randomized prospective trial versus traditional fiberglass cast

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OBJECTIVE — This study was designed to test the safety, effectiveness, and costs of off-loading with a novel, off-the-shelf removable device in the management of diabetic foot ulceration (DFU).

RESEARCH DESIGN AND METHODS — We prospectively evaluated off-loading of neuropathic plantar ulcers in 60 diabetic outpatients attending two diabetic foot clinics and compared healing rates at the 12-week follow-up, patient and severity of adverse events, healing time, costs and applicability of the device, and patients' satisfaction between those randomly assigned to total contact casting (TCC; group A) or to the Optima Diah walker (group B). Deep or infected ulcers were excluded.

RESULTS — No difference between groups A and B was observed in healing rates at 12 weeks (93 vs. 87%), healing time (6.5 ± 4.4 vs. 8.7 ± 3.4 weeks), and number of adverse events (six versus four). Treatment was significantly less expensive in group B, which showed a mean reduction of around 70% compared with group A ($P < 0.001$). Practicality was more favorable in group B, with a reduction of 77 and 20% of the time required for application and removal of the device, respectively ($P < 0.001$). Patients' satisfaction with the treatment was higher in group B ($P < 0.01$).

CONCLUSIONS — The Optima Diah walker is as safe and effective as TCC in the management of DFU, but its lower costs and better applicability may be of help in spreading the practice of off-loading among the centers that manage the diabetic foot.

Diabetes Care 30:380–390, 2007

Foot & Ankle International

AMERICAN ORTHOPAEDIC FOOT & ANKLE SOCIETY

Comparison of Removable and Irremovable Walking Boot to Total Contact Casting in Offloading the Neuropathic Diabetic Foot Ulceration

Alberto Piaggese, MD¹, Chiara Goretti, MD¹, Elisabetta Iacopi, MD¹, Giacomo Clerici, MD², Fabio Romagnoli, MD², Fabrizia Toscanella, MD⁴, and Cristiana Vermigli, MD³

Background: Despite its efficacy in healing neuropathic diabetic foot ulcers (DFUs), total contact cast (TCC) is often underused because of technical limitations and poor patient acceptance. We compared TCC to irremovable and removable commercially available walking boots for DFU offloading.

Methods: We prospectively studied 60 patients with DFUs, randomly assigned to 3 different offloading modalities: TCC (group A), walking boot rendered irremovable (i-RWD; group B), and removable walking boot (RWD; group C). Patients were followed up weekly for 90 days or up to complete re-epithelization; ulcer survival, healing time, and ulcer size reduction (USR) were considered for efficacy, whereas number of adverse events was considered for safety. Patients' acceptance and costs were also evaluated.


Results: Mean healing time in the 3 groups did not differ ($P = .5579$), and survival analysis showed no difference between the groups (logrank test $P = .8270$). USR from baseline to the end of follow-up was significant ($P < .01$) in all groups without differences between the groups. Seven patients in group A (35%), 2 in group B (10%), and 1 in group C (5%) (Fisher exact test $P = .0436$ group A vs group C) reported nonsevere adverse events. Patients' acceptance and costs were significantly better in group C ($P < .05$).

Conclusions: Our results suggest that a walking boot was as effective and safe as TCC in offloading the neuropathic DFUs, irrespective of removability. The better acceptability and lesser costs of a removable device may actually extend the possibilities of providing adequate offloading.

Level of Evidence: Level II, prospective comparative study.

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Diabetes Care

Diabetes Care N°30 03/2007, p. 586-590

An Off-the-Shelf Instant Total Contact Casting Device for the management of Diabetic Foot Ulceration 15/2016

FAI FOOT & ANCLE INTERNATIONAL

Comparison of Removable and Irremovable Walking Boot to Total Contact Casting in Offloading the Neuropathic Diabetic Foot Ulceration, Foot & Ankle International Vol 37, Issue 8, 2016

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